



Instructions For Use For

Oxevision Sleep

en-GB English

For use with Oxevision Sleep versions:

2.1 to 2.7

Oxevision Sleep System UDI-DI: 506075145SLEEP9J

Alternative electronic instructions for use are available at the website below, or you can request further paper copies from Oxehealth using the contact information at the end of these instructions.














Website: www.oxehealth.academy/instructions-for-use



Table of Contents

| | |
|--|-----------|
| Table of Contents | 2 |
| Symbols Used In This Manual And On Product Labels | 3 |
| Instructions For Use For Oxevision Sleep | 3 |
| Intended Use | 5 |
| Warnings | 5 |
| Cautions | 5 |
| Specification | 7 |
| Installation | 8 |
| Maintenance | 8 |
| How To Use Oxevision Sleep | 9 |
| The Overview Screen | 9 |
| Viewing the Sleep Report | 14 |
| Exporting a Sleep Report | 15 |
| Pausing Camera | 24 |
| Extending and Canceling a Pause | 24 |
| Changing Room Settings | 25 |
| Sending Feedback to Oxehealth | 27 |
| Manufacturer | 28 |

Symbols Used In This Manual And On Product Labels

| | | | |
|---|---|---|--|
|  | Read the Instructions for Use |  | Caution |
|  | Medical Device, according to the definition within Regulation (EU) 2017/745 on Medical Devices (as amended) |  | Unique Device Identifier |
|  | Lot or batch number (in the case of Vital Signs, the software version number) |  | Date of manufacture |
|  | Manufacturer |  | CE Marking, made in conformance with (MDR) EU 2017/745 for medical devices |
|  | UKCA Marking, made in conformance with the UK Medical Devices Regulations 2002 (as amended) |  | Authorised Representative in the EU. |
|  | Importer | | |

Instructions For Use For Oxevision Sleep



Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in inaccurate readings of sleep measurements.

These instructions for use are for the Oxevision Sleep software device, which is a product module within Oxehealth’s Oxevision Device. For instructions on how to use other product modules within the Oxevision system, refer to other instructions for use for those features. The Oxevision Sleep Device is intended to be used by healthcare professionals who have read this document.

R_x Only CAUTION: Federal law (U.S.) restricts this device to sale by or on the order of a licensed healthcare practitioner.



IMPORTANT - please read carefully

If any serious adverse incident occurs when using the Oxehealth Sleep medical device, where a subject or user is or could have been seriously harmed, it is vital that you inform Oxehealth using the contact details at the end of these instructions and you may also inform the FDA.

For details of how to inform the FDA of an adverse incident, please visit
<https://www.accessdata.fda.gov/scripts/medwatch/index.cfm>.

Intended Use

The Oxevision Sleep Device is an activity monitor designed and intended for documenting physical movements associated with applications in physiological monitoring.

The device's intended use is to analyze subject activity, movement and physiological sign data associated with movement during sleep and to extract information about certain sleep parameters from these movements and physiological sign data. The device provides a timeline of periods when a bed space is occupied, and periods when a subject is asleep when the bed space is occupied.

The Oxevision Sleep Device is software assessing video footage from a fixed-installation device for use within single occupancy bed spaces within hospitals, general care and secured environments.

The Oxevision Sleep Device is indicated for use on subjects 18 years of age or older.



Warnings

The Oxevision Sleep Device must not be used with any subject receiving critical care.

The Oxevision Sleep Device must not be used in rooms that contain or may contain:

- Domestic animals
- A shower or significantly sized water feature

Oxevision Sleep must not be used as the sole basis for making any decisions or recommendations relating to clinical diagnosis, subject care or treatment. Decisions or recommendations must be made by a clinician based on broader evidence than Oxevision Sleep alone provides.

It is possible that the introduction of the Oxevision Sleep Device could affect the behavior of individuals in your care. It involves the introduction of the Oxevision system including cameras and illuminators which could result in false perceptions of the capabilities, limitations and intended use of the software. It is essential that you engage with the individuals under your care and other interested parties to ensure that they fully understand the capabilities, limitations and the intended use of the software, using the training materials and subject information provided by Oxehealth. Failure to fully inform individuals could result in them becoming distressed.



Cautions

If the Oxevision Sleep Device fails to respond as described, discontinue use until the situation has been corrected. Contact Oxehealth for support through either the user interface, or the telephone number or address given at the end of this manual.

The Oxevision Sleep report may not display bed occupancy periods or sleep periods, or may display inaccurate bed occupancy or inaccurate sleep periods, in the following circumstances:

- Large airflows through the room.
- Other sources of movement in the room.

- If another living being is present in the room such as a carer or other person in or near the bed space (including hidden from view, for example by sheets).
- Human-like shapes are present in the bed space (e.g. a mannequin) or the bedding or other objects in the room resemble a human shape.
- Running, dripping or spraying liquids in the room.
- Moving patches of light or shadow in the image of the room as flashlights or lanterns.
- In the presence of sources of flashing or variable light, including other Near Infrared illumination, for instance from an existing CCTV installation.
- Sources of mechanical vibration in the image
- Vibration of the related hardware, caused for example by drilling into a nearby structure.
- Where insufficient body area is visible to the camera.
- The subject appears too bright in the camera image due to exceptionally bright illumination.
- The subject appears too dark in the camera image due to insufficient illumination, for instance if the infrared illuminators are obscured.
- The subject cannot be seen because the camera has been obscured.
- In subjects with pronounced sleep apnoea or who are deliberately controlling their breathing.
- The subject's torso is covered by a material with extremely low contrast edges.
- The subject's head is obscured (for example by a blanket or other bed clothes).
- The subject is in a dark area of the room, especially if a large part of the remainder of the room is covered in direct sunlight.
- The bed is moved from its original location in the room.
- In the 7-day period including the date of clock changes associated with Daylight Saving Time.

The Oxevision Sleep report may not display sleep periods, or may display inaccurate sleep periods, in the following circumstances:

- In subjects with low perfusion, anemia, poor circulation or blood volume deficiency.
- Where visible skin is significantly scarred or tattooed.
- Where visible skin has excessive make-up, sunscreen or other topical applications that hide the underlying skin applied.
- In subjects with tremor or motion disorders.
- In subjects under the influence of or withdrawal of alcohol.
- In subjects under the influence of or withdrawal of drugs.
- When subjects are talking (including vocalizations during sleep) or controlling their breathing.
- In the presence of sleep disorders, e.g. REM sleep disorders, periodic limb movement or parasomnias, or other disorders or syndromes leading to abnormal movement or breathing patterns during sleep.
- In the presence of cardiac arrhythmia or other forms of abnormal heart rhythm.

The Oxevision Sleep report displays data gathered from a specified bed space so may show data for any subject occupying that bed space.

Please be aware that:

- The device provides estimated sleep status. It is not a replacement for measurements of brain activity.
- The system is not intended to be used as an apnea monitor.
- The accuracy of the Oxevision Sleep System device has not been studied in patients with neuromuscular disorder.

The Oxevision Sleep device is not a real-time sleep monitor or a sleep apnea monitor. Its outputs do not indicate whether an individual is currently asleep, nor does it give any indication of the health or safety of an individual. The device cannot be used to determine whether an individual presents signs of abnormal or altered breathing.

The Oxevision Sleep device analyzes video captured from a video camera within a bed space illuminated by a Near Infrared Illuminator, and so obstruction of either of these parts of the device will severely compromise its operation. It is also possible that individuals will attack or attempt to damage the housing containing the camera and illuminator, and this can also result in compromising the operation of the system and the software. Ensure that the housing containing the camera and illuminator is kept free from obstruction according to the maintenance directions in the Oxevision Sleep Instructions for Use, and inform Oxehealth of any attack on or interference with the housing.

Oxehealth software and related hardware installed inside or outside the room environment is not to be installed or maintained by anyone other than an Oxehealth trained representative. Any attempt by a non-Oxehealth representative to interfere with the hardware, operating system administration or software configuration outside that described in this manual may result in damage to the device performance.

If you make significant physical changes to a room, for instance visible structural changes or adding large items of furniture, this may affect the device performance and you should inform Oxehealth.

Specification

| Clinical-related evaluation of Oxevision Sleep Device outputs | | | |
|---|---------------------|--------------------|--------------------|
| Bed State outputs | Sleep State outputs | | |
| Accuracy | Overall Agreement | Positive Agreement | Negative Agreement |
| >95% | >82% | >88% | >55% |

Installation

Installation will be carried out by an Oxehealth approved installer, according to installation instructions and verification procedures provided by Oxehealth.

Maintenance

The Oxevision Sleep device analyzes video from an infrared illuminated room that has been captured by a video camera to measure activity, movement, the pulse rate and breathing rate (chest wall movements) of a subject. Hardware is often installed within a housing and so the cleanliness of the housing must be maintained. When cleaning the housing unit and/or touch screen devices, do not use solvent based, ammonia, caustic or abrasive cleaners or abrasive cleaning materials on the clear front of the housing unit or on touch screen device displays, as this may scratch or fog the material, and potentially degrade the ability of the device or the user's ability to see the subject clearly and may reduce the reliability and accuracy of the readings.

If components within the sealed housing are directly exposed to liquids, vapors or dust then they may be damaged and cease to operate correctly. If this happens, stop using the Oxevision Sleep software device and inform Oxehealth.

In the event of failure of any part of the Oxehealth product, inform Oxehealth.

Only Oxehealth or Oxehealth trained maintenance personnel may make modifications to this device as these may affect device performance. If you are in doubt whether changes could affect the device, you should contact Oxehealth who can advise.



IMPORTANT SAFETY AND SECURITY NOTICE - please read carefully

Maintaining the privacy of subjects, and protecting the confidentiality of their data is of the utmost importance. Ensure that you observe all of your organization's safety, privacy and security policies and guidelines, and do not relocate any equipment that has been installed for the Oxehealth software, or change any security or access arrangements to the equipment without first informing Oxehealth.

If you are using any mobile devices to access your Oxehealth software, pay particular attention to ensuring these devices are used only by authorized and trained staff, using the devices as intended.

Oxevision Sleep is not intended to be used as a surveillance tool, and you must use the functions to view the video feed in accordance with your organization's privacy and security policies. Failure to do so puts you and your organization at risk of violating individuals' rights, and data privacy laws.

NOTE

If you are using any mobile devices to access your Oxehealth software, they should be kept charged. In the event that these devices cannot be accessed or have a depleted battery, Oxevision Sleep can be accessed using the fixed installation screen.


How To Use Oxevision Sleep

NOTE

In the case for this Instructions For Use, "Room 2", "Room 10", "Bedroom 4", etc., are customizable names for the areas in which Oxevision is installed. References to "rooms" in the text can also refer to the locations, or bed spaces where Oxevision is installed.

The Overview Screen

The **Overview Screen** shows all of the rooms connected to the Oxevision Sleep device:



IMPORTANT - please read carefully

Depending on the number of rooms and layout chosen, you may need to scroll down to see all of the room cards.

Each room is shown as a clickable panel on the Overview Screen, by clicking a room, this brings up the **Room Dashboard** for that room:

The screenshot displays the Oxevision Room Dashboard for Room 5. At the top, there is a navigation bar with the Oxevision logo, a menu icon, and the text 'Observation Round' and '12 Dec 2024, 10:52'. Below the navigation bar, the breadcrumb 'Home > Room 5' is visible. A note states: 'Note: The below previews do not include all report information. For full information, click on view report.' The main content area is divided into several sections:


- Activity Report:** A bar chart showing activity levels from 22:00 to 10:00. The chart shows a green bar from 22:00 to 00:00, indicating activity during the night. A 'View Report' link is present.
- Alerts:** A table listing alerts with columns for Time and Type. The alerts are:

| Time | Type |
|-------|-----------------|
| 10:20 | In Bathroom |
| 10:14 | Multiple People |
| 05:51 | In Bathroom |
| 21:14 | In Bathroom |
- Summary:** A table showing activity and average values:

| Activity | Average |
|-------------|----------|
| In Bed | 11h 4min |
| In Bathroom | 1x |
| Room Entry | - |
- Room Information:** A card showing 'Room In Bed 1 hr 6 mins 5 secs' and 'Occupied by: Hank Seymour'. It includes buttons for 'Pause Camera for 15 mins' and 'Patient Configuration'.
- Vital Signs:** A card with a 'Take Vital Signs' button. It shows 'Last vital signs recorded: 14 hrs 49 mins 37 secs' and buttons for 'Pulse Rate' and 'Breathing Rate'.

On the right side of the dashboard, there is a vertical sidebar with 'Alerts 2' and 'Feedback' buttons.

Note that you may need to scroll down on the **Room Dashboard** to view the Sleep report:



Note that the room **Room Dashboard** includes the **Sleep Report**, and the Oxehealth Vital Signs features. Depending on your software configuration of the room, your view may vary and the **Room Dashboard** may include other menu items for other product modules, for example Viewing Vital Signs Trends (refer to other Oxehealth Oxevision product module instructions for use for information).


On the **Overview Screen**, the large icons indicate the status of the room as follows:

- **Connection Lost:** the Oxevision Sleep device is not connected to the room.




The main text in the room icon will display the nature of any software disconnection, and the length of time that the connection has been down.


In this state, it will not be possible to access the **Sleep Report** or other actions on the Room Dashboard.



- **Connected:** the Oxevision Sleep device is connected to the room and running, such that the software is receiving video and making estimations.



The main text in the room icon will display information from other modules about what is happening in the room, and the time that the room has been in this state. Other modules may also modify the color of the room icon while they are processing data, or when they determine the room to be empty.




Refer to other product module instructions for use for further information. In this state, it will be possible to access **Sleep Reports**.

IMPORTANT - please read carefully

When other Oxehealth product modules are running, additional information about what is happening in the room might be shown. The Oxevision Sleep is still operational, and it is still possible to access Sleep Reports.


- **Camera Off:** the Oxevision Sleep device is connected to the room, but a user has chosen to stop measuring, such that the Oxevision Sleep device has ceased receiving video and making estimations for this room.



In this state, it will be possible to access the **Sleep Report**. In this state, it will not be possible to use the **Pause Camera** function on the **Room Dashboard**.

- **Camera Paused:** the Oxevision Sleep device is connected to the room, but a user has chosen to pause the camera measuring, such that the Oxevision Sleep device has ceased capturing video and making estimations temporarily for this room.


The time for the camera pause will continuously update, and is a countdown timer, indicating how long is left before Oxevision Sleep automatically restarts capturing video and processing.



In this state, it will be possible to access the **Sleep Report** and use the **Extend Pause**, **Cancel Pause** **Camera**.

Viewing the Sleep Report

To view your sleep report for the bed space, click the room icon for the room the subject is in, then select the **View Report** button on the **Sleep Report** from the **Room Dashboard**. You will then see the sleep report screen below:



The Sleep Report shown above includes 24-hour timelines showing a subject’s sleep and bed occupancy patterns, and the operational status of the system over a reporting period (i.e. last seven days). Each timeline will include the times when a subject was detected to be in the bed space, and the times a subject was asleep for those periods when a subject was detected in the bed. Times at

which the system has been switched off, or when bed occupancy or sleep status could not be computed are also presented.

On the right hand side of the screen are a set of key statistics about the current reporting period. It will include the number of hours the subject was in-bed and the number of hours the subject was asleep for each 24-hour timeline.

The time of the day at which a Sleep Report is updated with new data may vary depending on your software configuration of Oxevision Sleep.

To go back to the **Room Dashboard** page, click the Room button at the top left hand corner of the screen.

To return to the homescreen, click the **Home** button at the bottom left hand corner of the screen.

Exporting a Sleep Report

To export a Sleep Report you are viewing for the given room, click the **Export** button at the top right of the screen:



You will then see the **Choose Recipient** screen below:

A screenshot of a web application interface. At the top left is the "oxevision" logo and a menu icon. At the top right is the date and time "19 Sep 2024, 12:29". Below the header is a breadcrumb trail: "Home > Room 4 > Sleep Report > Choose Recipient". The main content area is divided into two columns. The left column has an "Email" section with a text input field containing "Recipient's email" and "@oxehealth.com". Below that is a "Reason For Export" section with a note "Do not include personally identifiable information, e.g. names, dates of birth, etc." and a large text area with the placeholder "Type your reason here...". At the bottom of this section, it says "Characters remaining: 255" and "This email will also be sent to the following default recipients for audit purposes: some.email@oxehealth.com". The right column has a "Summary" section with three sub-sections: "Exporting" (Sleep Report), "For" (Room 4), and "Time Period" (12 Sep 2024 - 19 Sep 2024). On the far right is a vertical sidebar with "Alerts 2" and "Feedback" buttons. At the bottom right of the main content area is a "Confirm Export" button.

Enter the email address you wish to send the report to in the **Email** field, and describe a reason for exporting in the **Reason for Export** field. These fields must be completed for the **Confirm Export** button to become available and to successfully send the report.

You will only be able to send the report to a recipient with an email address within your organization, and you only need to enter the part of the email address before the '@' symbol. When you export a report, a copy may also be sent to other recipients that your organization has chosen for audit purposes. The email address of this person will be indicated on the **Choose Recipient** screen shown above.

NOTE


Do not include any subject identifying information in the **Reason For Export** section, e.g. names, date of birth.

NOTE


Please take care in entering email addresses, as Sleep Reports may contain healthcare data. Please ensure that this data is shared according to your organization's policies.

To cancel this action from the **Choose Recipient** screen, click the **Sleep Report** button at the top of the screen to return to the Sleep Report screen or the **Home** button to return to the **Overview screen**.

When you are ready to export the report, click the **Confirm Export** button. A message will appear which reads "Submitted. A report email will be with you in the next 5 minutes" as shown below:



If you are operating a tablet device and this has lost connection to WiFi, or the system encounters other errors attempting to send the email, you will see an error message as shown below:




To ensure that Sleep Reports can be sent, move to a place where your tablet device has a good Wifi connection and click **Confirm Export** again. If the problem persists, please contact Oxehealth using the details on the screen.

Once successful, an email will be sent to the address you provided, with the report attached as a PDF file. An additional copy of this email will also be sent to the other email addresses shown on the screen for audit purposes. The Sleep Report will contain one page which will display one row of sleep data per day for the selected room.


Emails should be received by the recipients, and default recipient(s) for the device group, within 5 minutes.

If you have not completed the **Email** and **Reason for Export** fields, error messages will be shown and the report will not be sent. After successfully exporting, you will then be returned to the Overview Screen.

Alternatively, you can export Sleep Reports for a number of rooms over a selected date range by clicking the **Sleep Report** link in the navigation menu on the **Overview Screen**:



You will then see the **Select Rooms** screen shown below, allowing you to choose multiple rooms to export the reports for:



All rooms available to the device will be shown as multi-select options within the **Select Room** page. One or more rooms can be selected by clicking on the room name button. Clicking the **Select All** checkbox will select all rooms available. Selected rooms will be displayed in a solid (green) color. Unchecking the **Select All** checkbox will deselect all rooms.

The **Choose Dates** button will be disabled until at least one room is selected.

To go back to the **Exports** screen to change the reports being exported, click the **Exports** button at the top left.

After selecting the rooms you would like to export Activity Reports for, click **Choose Dates**, and you will then see the **Choose Dates** screen shown below:

oxevision® 19 Sep 2024, 12:34

Home > Select Rooms > Choose Dates

You are exporting Sleep Report for 2 Rooms

Start End

Home Choose Recipient

Alerts 2 Feedback

The time frame for the export can be selected using the **Start** and **End** date dropdown boxes from the **Choose Dates** page.

The **Choose Recipient** button will be disabled until at least one date is selected.

The maximum date range is 14 days. If invalid dates are selected, an error message will be displayed and the **Choose Recipient** button will continue to be inactive.

To go back to the **Select Rooms** screen to change the rooms being exported, click the **Select Rooms** button at the top left.

Once you have selected the date range you would like to export, click the **Choose Recipient** button and you will then see the **Choose Recipient** screen below:

Enter the email address you wish to send the reports to in the **Email** field, and describe a reason for exporting in the **Reason for Export** field. These fields must be completed for the **Confirm Export** button to successfully send the reports.

You will only be able to send the report to a recipient with an email address within your organization, and you only need to enter the part of the email address before the '@' symbol. When you export a report, a copy may also be sent to other recipients that your organization has chosen for audit purposes. The email address of this person will be indicated on the **Choose Recipient** screen shown above.

NOTE

Do not include any subject identifying information in the **Reason For Export** section, e.g. names, date of births.


To go back to the **Choose Dates** screen to change the dates being exported, click the **Choose Dates** button at the top left. To cancel this action from the **Export Report** screen, click the **Home** button to return to the **Overview** screen.

A summary of the export is displayed on the right of the **Choose Recipient** screen which shows details of:


- **Exporting** - The type of report being exported;
- **For** - The rooms selected and;
- **Time period** - The time frame selected.

The export options can be edited by selecting the edit icons within the Summary section. Selecting an edit icon will return you to the user to the respective **Exports**, **Select Rooms** or **Choose Dates** screen. Follow the earlier instructions for editing the information on these screens.

When you are ready to export the report, click the **Confirm Export** button. A message will appear which reads “Submitted. A report email will be with you in the next 5 minutes”. as shown below:



If you are operating a mobile device and this has lost connection to WiFi, or the system encounters other errors attempting to send the email, you will see an error message as shown below:



To ensure the reports can be sent, move to a place where your mobile device has a good Wifi connection and click Confirm Export again. If the problem persists, please contact Oxehealth using the details on the screen.

Once successful, an email will be sent to the address you provided, with the reports attached as individual PDF files (one per room). An additional copy of this email will also be sent to the other email addresses shown on the screen for audit purposes. Each Sleep Report will contain up to 7 days of data per page.

Emails will be received by the recipient, and default recipient(s) for the device group, within 5 minutes.

If you have not completed the **Email** and **Reason for Export** fields, error messages will be shown and the report will not be sent. After successfully exporting, you will then be returned to the **Overview Screen**.

Pausing Camera

To pause any room, click the room icon for the room the subject is in, then select **Pause Camera for 5 Minutes** from the **Room Dashboard**.

NOTE

The duration of the camera pause may not be 5 minutes if the default pause camera duration for this room has been changed. For instance, if the default pause camera duration has been changed to 45 minutes, the button will read **Pause Camera for 45 Minutes**.

Extending and Canceling a Pause

If you want to keep a room paused, or cancel the action and reconnect the room, click the room icon for the room the subject is in, and the **Room Dashboard** will look as it does below:

The screenshot displays the Oxevision interface for Room 5. At the top, it shows 'Observation Round' and the date '12 Dec 2024, 11:04'. The main content area is divided into several sections:

- Activity Report:** A bar chart showing activity levels from 23:00 to 11:00. A green bar indicates activity from 23:00 to approximately 07:00.
- Alerts Table:**

| Time | Type | Action |
|-------|-----------------|--------|
| 10:20 | In Bathroom | ▶ |
| 10:14 | Multiple People | ▶ |
| 05:51 | In Bathroom | ▶ |
| 21:14 | In Bathroom | ▶ |
- Summary Table:**

| Activity | Average |
|-------------|----------|
| In Bed | 11h 4min |
| In Bathroom | 1x |
| Room Entry | - |
- Room Status:** 'Camera Paused 9 mins 19 secs'. Occupied by: Hank Seymour.
- Buttons:** 'Cancel Pause Camera', 'Extend Pause 15 mins', and 'Patient Configuration'.
- Vital Signs:** 'Take Vital Signs' button. Last vital signs recorded: 15 hrs 1 min 7 secs.

To cancel the pause and reconnect the room, click **Cancel Pause Camera**.

To extend the pause for another 15 minutes, click **Extend Pause 15 Minutes**. The camera will now remain paused for the time that previously remained, plus the time added.


Click **Home** to return to the **Overview screen**


NOTE

The extension of the pause may not be 15 minutes if the default pause camera duration for this room has been changed. For instance, if the default pause camera duration has been changed to 45 minutes, the button will read **Extend Pause 45 Minutes**.

Changing Room Settings



To enable or disable Oxevision Sleep device software processing in any of the rooms, or change the default pause camera duration, click the **Settings** link from the navigation menu.:





On the settings screen, you can control which rooms are connected, and set the default pause duration.

This screen shows the on or off status for each room, which is operated as follows:


| | |
|---|---|
| <p>Room 1 </p> | <p>The room is chosen to be on, such that the software is running, receiving and processing video.</p> <p>To turn the room off (stop the software from running, such that the software is not receiving or processing video), click the green toggle switch. The toggle switch should then look like the picture below.</p> |
| <p>Room 1 </p> | <p>The room is chosen to be off, such that the software is not running, receiving and processing video.</p> <p>To turn the room on (start the software, such that the software is receiving and processing video), click the black toggle switch. The toggle switch should then look like the picture above.</p> |

To change the **Default Pause Camera** duration, click the up arrow above the pause camera duration to increase the pause camera duration **minutes** and click the down arrow below the pause camera duration to decrease the pause camera duration **minutes**.

Click **Home** to return to the **Overview screen**.

Sending Feedback to Oxehealth

To send feedback of any sort to Oxehealth, click the **Feedback** tab in the top right corner of the screen. You will then see the **Feedback screen** below:



Select the **Feedback Type** (if this is available) and fill in the **Name**, **Contact Number** and **Message** text field with the information you would like to share, your issues, suggestions or questions and click **Send** to send the information to Oxehealth. Please provide your name and contact details,, as this will help a support representative to contact you if you require assistance. We will then receive your message automatically, and you will be returned to the **Overview screen**.

If you do not wish to send feedback, click on the **Feedback tab** and the feedback form will close



CAUTION

Do not include any subject identifying information in the feedback, e.g. names, date of birth.

Manufacturer



Oxehealth Limited

Bee House
Eastern Avenue
Milton
Abingdon
England
OX14 4SB

Telephone: US: +1 (888) 871-8108

Telephone UK: 0800 030 6781

Website: www.oxehealth.com

Email: support@oxehealth.com